

# LABCONTROL



**EASYPAB THE NEW GENERATION**  
**Controlling the ventilation in lab environments**

**TROX<sup>®</sup> TECHNIK**  
The art of handling air

Paul Absolon

Trox Belgium

Labcontrol Product Manager

**DIN 1946 Part 7 Lab ventilation**

**DIN EN 14175 Fume hoods**

**Different Lab Ventilation Layouts**

**Requirements for Labcontrol Systems**

## DIN 1946 Part 7 Lab ventilation

minimum exhaust according to DIN 1946,  
P7 ( $25 \text{ m}^3/\text{h}/\text{m}^2$ )

### **DIN EN 14175 Fume hoods**



**No Fixed Value is given for :**

**Min Volume Flow**

**Min Face Velocity**


**Manufacturer Guaranteed**

## Fume Hoods compensated by Supply Air



x = Konstante Differenz Zuluft/Abluft  
zur Sicherstellung der Überströmung

## Fume Hoods comp Supply/Extract Air




$x = \dot{V}$ -Differenz zur Einhaltung des Unterdrucks

### Which requirements can be fulfilled with LABCONTROL Systems?

- **Fast Control Circuits**
- **Volume Flow Control**
  - With temperature control concepts
  - Room balance concepts with many Master units
- **Room Pressure Control**
  - Volume flow pressure cascade with:  
Duct pressure control / Room pressure control
- **Fume Hood Control**
  - variable with distance sensor, 2- or 3-point control,  
Face velocity sensor
  - constant volume flow
- **Monitoring concepts**
- **Easy installation and wiring**
- **Far distance maintenance**

Next generation of the Controller hard- and software



**EASYLAB - TCU3**



This wishes should be fulfilled

- **EASY to adapt**
- **EASY to install**
- **EASY to design**
- **EASY to use**
- **EASY to commission**

## The modulare priciples make it possible:

- Volume flow controller with baffle measurement or according to the venturi principle
- Plug and Play connection between the controllers
- Adaptable User terminals
- Flexible Adaptation of the digital communication
- Controller connections easy to reach
- Possible UPS
- Failure diagnostics very easy to handle
- Flexible power supply 24V AC / 230V AC
- If needed auto zero function possible

*Lets have a closer look at those points ...*


## ● EASY to install



### Terminal assignment out side:

- 1 LED for failure indication
- 2 Connection for sash window contact EN 14175 (500mm)
- 3 Connection user terminal 1
- 4 Connection user terminal 1
- 5 Connection Actuator
- 6 Connection face velocity sensor (VS-TRD)
- 7 Communication line
- 8 Communication line
- 9 Cable bushing on both sides
- 10 Connection fume hood light (expansions modul EM-LIGHT)

- Max. 24 controllers in a EASYLAB-System
- Controller connection over a communication line KL
- Only line network topology max. 300 m
- Termination on the start / end of line necessary!



- Plug-In Standard-Patchkabel (green)  
Cable typ: S-FTP

## EASYPAB Regler TCU3

**A Fast Runner is needed**

**Precision Motor (Belimo)**

**0-90° <= 3 s**

**Accuracy 0,2°**

**8 Nm**

**Is use for all elab vav systems**

**Exception :**

**TVJ/TVT big vav systems (15Nm)**




### Overview

- **EM-TRF / EM-TRF-USV**  
Power supply 230 V / USV function
- **EM-LON**  
LonWorks<sup>®</sup>-Gateway
- **EM-AUTOZERO**  
automatic zero point adjustment
- **EM-LIGHT**  
connection for fume hood light



## EASYLAB expansion module

- **EM-TRF** power supply 230 V AC
- **EM-TRF-USV** power supply 230 V AC  
+ USV function incl. Accu
- **Mounting in the factory or retrofit**



## EM-LON

- **LonWorks<sup>®</sup>- gateway FT10**  
Free Topology Network
- **Automatic detection of the module**  
no configuration necessary!
- **Mounting in the factory or retrofit**






### EM-AUTOZERO

automatic zero point adjustment


- **Mounting in the factory or retrofit**
- **Plug-In electrical connection**
- **Automatic detection of the module**  
no configuration necessary!



## EM-LIGHT

connect for fume hood light

- For fume hood controller mounting in the factory or retrofit
- Matching plug not in scope of delivery






# Adapting user terminals for fume hood

**A lot of additional!**

- Set and current values
- Alarm and fault indication

**EASY to use**




**New terminals – the next generation**

**Complete new room scenes can be done**




## The room terminal:

A lot of extra information and modes can be modified!



Can be easy installed:



## Take you right to work whenever you want!




**Complete new solutions:**

**Hand-mode**

**Enables the possibility of from  
the BMS independent user modes  
(if wanted with timer!)**


*That is fast!!*

**The wiring can be done within minutes!**



**EASYLAB compatible**

**EASYLAB can for sure be combined with all volume flow controller series!**





***Complete new development !!!***

**TVLK –  
high precision injection moulded part with baffle**



*Complete new development !!!*

**TVLK – and as venturi option**



***Complete new development !!!***


**TVLK – venturi option**

**Why should you miss the cleaning possibility?**



# Everything started with the fume hood ...



## The fume hood – the start of LABCONTROL



**A source of possibilities**



For sure we can do this ...

## Fume Hood with two volume flows





... this for sure as well ...

## Fume hood with 3 volume flows switchable




... often copied but never reached ...

## Face velocity control



Inklusive: ...


# LABCONTROL VS-TRD



**Circuit board outside  
the airstream**


**Measuring unit**

**room air inlet from  
downside =>  
better resistant against  
room turbulences**






**... Reaction on thermal loads!**



**Adapted voöume flow to Vmax**


**Heat source**



**Alternative:  
Smoke or  
fire detection**



... If this is wanted we can do it as well!

## Distance sensor – We can do as well!





Simply have a look on the costs!

The linear function ...



Simply a little bit safer!  
The optimized function ...





➤ Some screenshots for fume hood config

## Device - Basic configuration

Device name	TroxBelgium	DemoRoom	Fumehood	ID -	1
C-value	25.4	based on l/s	You will find the valid C-value for this volume flow controller on the device product label.		
Technical Vmin	198	m <sup>3</sup> /h	(lower volume flow limit of volume flow controller)		
Technical Vnominal	1296	m <sup>3</sup> /h	(upper volume flow limit of volume flow controller)		

Geräteserie: TVLK  
Gerätegröße: 250-D16  
Gerätefunktion : FH-VS

Vnenn: 702 m <sup>3</sup> /h	Ji
Pw: 215 Pa/Vnom	2
<b>(13,3)</b>	S
Vmin: 108 m <sup>3</sup> /h	

## Load configuration settings from file

Select a file to load existing configuration settings.

Load configuration file

Continue to transmit configuration settings

## Operation mode preset for this fume cupboard controller

Room operation mode presets will be taken over

### Modification of room operation mode

- Room mode High is changed to Fume cupboard Standard mode
- Room mode Shut off is changed to Fume cupboard Standard mode
- Room mode Standard is changed to Fume cupboard Shut off mode

- Room mode Low is changed to Fume cupboard Standard mode
- Room mode Low is changed to Fume cupboard Shut off mode

Room operation mode presets will be ignored / Local presets are possible

## Control concept for fume cupboard control

- FH-VS Face velocity sensor at AI5  
Control concept to guarantee a defined face velocity
- FH-DS Sash distance sensor at AI4  
Linear control concept
- FH-DV Sash distance sensor at AI4  
Security optimized control concept
- FH-3P 3 point control via 2 switch contacts at DI2 and DI3
- FH-2P 2 point control via one switch contact at DI2
- FH-F Fixed value control



## Equipment components


- Alarm Sash monitoring EN 14175 Via switching input DI1
- Fume scrubber (DI4)
- Supportive flow technology (DI5)
- Automatic sash device (Third-party supplier)  
(While using TROX FSE this configuration should be deactivated)
- Motion detector
  
- Smoke extraction function
  
- Operation mode dependent relay switching (e.g. fans)
  
- Integration of switchable constant volume flows into room balance  
(Hoods, snorkels, etc. with switching output)
- Integration of variable volume flows into room balance  
(Hoods, snorkels, etc. with 0-10V analogue output)

Cancel

Previous

Next

## Integration of switchable constant volume flows into room balance



Constant volume flow (D11)

Constant volume flow (D12)

101 m<sup>3</sup>/h

Exhaust air

Make contact = Integrates constant volume flow

Constant volume flow (D13)


Constant volume flow (D14)

Constant volume flow (D15)

Constant volume flow (D16)



## Integration of variable volume flows into room balance



Variable volume flow (A12)

Exhaust air ▾

Characteristic

min 2.00 V

equals

min 0 m³/h

max 10.00 V

equals

max 202 m³/h

Signal smoothing (1-50)

50

Variable volume flow (A13)

Supply air ▾

Characteristic

min 2.00 V

equals

min 0 m³/h

max 10.00 V

equals

max 299 m³/h

Signal smoothing (1-50)

50

Variable volume flow (A14)

Variable volume flow (A15)

## Control panel



Optical alarm

Red alarm - continuous



Acoustic alarm

Alarm duration limited to

15

s



Sash monitoring selected

Acoustic alarm duration limited to

Enable key - High mode

High mode unlimited



Enable key - Low mode



Enable key - Shut off mode



Automatic sash device not selected



Enable key - Fume cupboard light

Switch off fume cupboard light when activating Low Mode or Shut off mode



Enable key - Hand mode

## Control panel - Display options

Enable display options of BE-LCD

Enable display options of BE-SEG / BE-LCD

Display language

English

Display setpoint/current volume flow values

Display current face velocity value

Volume flow unit


m<sup>3</sup>/h

**Control concept: Face velocity transducer - Control strategy to guarantee a defined face velocity**

**Volume flow set values**

Standard mode / Vmin	<input type="text" value="202"/>	m³/h
Standard mode / Vmax	<input type="text" value="698"/>	m³/h
High mode	<input type="text" value="698"/>	m³/h
Low mode	<input type="text" value="202"/>	m³/h

Consider diversity control (Limit volume flow to allowed maximum of Total room exhaust air)



**Face velocity**

Setpoint	<input type="text" value="0.50"/>	m/s								
Tolerance around setpoint (+/-)	<input type="text" value="0.05"/>	m/s	Characteristic	min	<input type="text" value="2.00"/>	V	equals	min	<input type="text" value="0.00"/>	m/s
				max	<input type="text" value="10.00"/>	V	equals	max	<input type="text" value="1.00"/>	m/s
			Signal smoothing (1-50)		<input type="text" value="50"/>					
Time target - Setpoint adaption open sash	<input type="text" value="3"/>	s								
Time target - Setpoint adaption close sash	<input type="text" value="7"/>	s								

## Controller alarm settings

Category 1 (Smoke extraction / UPS Battery operation)

Activate monitoring

Category 2 (Volume flow control functions)

Volume flow monitoring (Standard mode, High mode)

Volume flow monitoring (Low mode)

Alarm delay time  s

Face velocity monitoring (Standard mode, High mode)

Trigger alarm below  m/s

Category 3 (Hardware failure)

Activate monitoring

## Service interval

Specified service period  days

Period since last service  days

Signaling at end of service interval

## Volume flow measuring point

- Use internal volume flow transducer
- Use external volume flow transducer at terminal A11 (e.g. Ex-Controller)

Characteristic min  Pa  
max  Pa

Signal smoothing (1-50)

Attention: Factory settings  
Values are dependent on used transducer.

Modification only by instructed technical staff !

## Volume flow control

Control tolerance  % with additional absolute tolerance of  l/s equals 14 m<sup>3</sup>/h

- Automatic calculation of limit value to achieve maximum actuating value

Open damper 132 m<sup>3</sup>/h  
Close damper 132 m<sup>3</sup>/h

Open damper

Minimum actuating value  [0.01 ... 100] Maximum actuating value  [0.01 ... 100]

- Manual adaption of limit to achieve maximum actuating value.

Close damper

Minimum actuating value  [0.01 ... 100] Maximum actuating value  [0.01 ... 100]

# Fume Hood Config

## Analogue output A01 - Current volume flow (Fume cupboard)

Characteristic	min	<input type="text" value="2.00"/>	V	equals	min	<input type="text" value="0"/>	m <sup>3</sup> /h
	max	<input type="text" value="10.00"/>	V	equals	max	<input type="text" value="1296"/>	m <sup>3</sup> /h

## Analogue output A02 - Total volume flows (Room)

- Total exhaust air      Fume cupboards / Room exhaust air / Integrated exhaust air / Constant exhaust air
- Total supply air      Room supply air / Integrated supply air / Constant supply air
- Setpoint supply air      Total exhaust air minus integrated and constant supply air and consideration of a volume flow balance

Constant exhaust air, constant supply air as well as volume flow balance values require a configuration within a controller with Room-Management-Function.

Characteristic	min	<input type="text" value="2.00"/>	V	equals	min	<input type="text" value="0"/>	m <sup>3</sup> /h
	max	<input type="text" value="10.00"/>	V	equals	max	<input type="text" value="0"/>	m <sup>3</sup> /h

## Analogue output A03 - Damper position

Characteristic	min	<input type="text" value="2.00"/>	V	equals	0 % Damper position (Damper closed)
	max	<input type="text" value="10.00"/>	V	equals	100 % Damper position (Damper open)

## Analogue output A04 - Actuator

Characteristic	min	<input type="text" value="2.00"/>	V	equals	0 % Damper position (Damper closed)
	max	<input type="text" value="10.00"/>	V	equals	100 % Damper position (Damper open)



## Room controller

- Integration of switchable constant volume flows into room balance (Hoods, snorkels, etc. with switching output)
- Integration of variable volume flows into room balance (Hoods, snorkels, etc. with 0-10V analogue output)
- Diffusor volume flow optimisation

## Room-Management-Function (RMF)

- Activate Room-Management-Function (RMF) (Centralised room configuration and integration for room operation mode presets, shift signals, etc.)
- Support EASYLAB Room control panel
- Sun blinder control
- Room operation mode preset via switching inputs (DI)
- Shut off for all room exhaust controller within Low Mode
- Diversity Monitoring / Diversity control
- External volume flow shift Volume flow shift via LonWorks<sup>®</sup>/BACnet/Modbus ▼
- Pressure control internal ▼  
with two alternative setpoints (switchable via LonWorks<sup>®</sup>/BACnet/Modbus) ▼
- Use door contact DI5
- Room exhaust air optimisation



### Volume flow split-up if several room controllers

Within RMF selected mode: Automatic split-up 100 %


Only if manual volume flow split-up is selected within RMF:

Supply air fraction of this controller  %

### UPS function

Mode during supply voltage failure

- Keep up normal controller action
- Special operation mode - Open mode
- Special operation mode - Shut off mode
- Keep operation mode / Keep damper position



### Controller alarm settings

Category 1 (UPS Battery operation)  Activate monitoring

Category 2 (Volume flow control functions)  Activate volume flow monitoring (Standard mode, High mode)

Activate volume flow monitoring (Low mode)

Alarm delay time  s

Category 3 (Hardware failure)  Activate monitoring

### Diffusor volume flow optimisation

Shut off 1 (DO5) if volume flow is below volume flow limit 1  m<sup>3</sup>/h      Switching tolerance 1  %

Shut off 2 (DO6) if volume flow is below volume flow limit 2  m<sup>3</sup>/h      Switching tolerance 2  %

## Analogue output A01 - Current volume flow (controller)

Characteristic	min	<input type="text" value="2.00"/>	V	equals	min	<input type="text" value="0"/>	m³/h
	max	<input type="text" value="10.00"/>	V	equals	max	<input type="text" value="1458"/>	m³/h

## Analogue output A02 - Total volume flows (Room)

- Total exhaust air      Fume cupboards / Room exhaust air / Integrated exhaust air / Constant exhaust air
- Total supply air      Room supply air / Integrated supply air / Constant supply air
- Setpoint supply air      Total exhaust air minus integrated and constant supply air and consideration of a volume flow balance

Constant exhaust air, constant supply air as well as volume flow balance values require a configuration within a controller with Room-Management-Function.

Characteristic	min	<input type="text" value="2.00"/>	V	equals	min	<input type="text" value="0"/>	m³/h
	max	<input type="text" value="10.00"/>	V	equals	max	<input type="text" value="0"/>	m³/h

## Analogue output A03 - Damper position

Characteristic	min	<input type="text" value="2.00"/>	V	equals	0 % Damper position (Damper closed)		
	max	<input type="text" value="10.00"/>	V	equals	100 % Damper position (Damper open)		

## Analogue output A04 - Actuator

Characteristic	min	<input type="text" value="2.00"/>	V	equals	0 % Damper position (Damper closed)		
	max	<input type="text" value="10.00"/>	V	equals	100 % Damper position (Damper open)		

## Room operation mode preset via switching contacts (Digital inputs DI1-DI6)

- |  |               |                                   |
|--|---------------|-----------------------------------|
| <input checked="" type="checkbox"/> Room operation mode preset via DI1 | Shut off mode | Make contact = Activates function |
| <input checked="" type="checkbox"/> Room operation mode preset via DI2 | Low mode      | Make contact = Activates function |
| <input checked="" type="checkbox"/> Room operation mode preset via DI3 | Standard mode | Make contact = Activates function |
| <input type="checkbox"/> Room operation mode preset via DI4            |               |                                   |
| <input type="checkbox"/> Room operation mode preset via DI5            |               |                                   |
| <input type="checkbox"/> Room operation mode preset via DI6            |               |                                   |
- If no valid room operation mode preset is given by switching contacts, the preset is Standard mode.

## Volume flow settings for the room

### Exhaust air

Minimum Total room exhaust air

High mode

m<sup>3</sup>/h

Standard mode

m<sup>3</sup>/h

(compliant with DIN 1946 Part 7 25m<sup>3</sup>/h per m<sup>2</sup> Lab floor space)

Low mode

m<sup>3</sup>/h

Total amount unmeasured Constant exhaust air

m<sup>3</sup>/h

(e.g. Exhaust by RN controller)

### Supply air

Total amount unmeasured Constant supply air

m<sup>3</sup>/h

(e.g. Supply air by RN controller)

### Volume flow split-up

Automatic split-up for Exhaust air volume flows

(to all exhaust controllers of room)

Automatic split-up for supply air volume flows

(to all supply air controllers of room)


## Diversity Monitoring / Diversity control

Total room exhaust air limit  m<sup>3</sup>/h  
(defined by air conditioning plant layout)

- Monitor the adjusted limit
- Control the adjusted limit (Selective diversity control)

Tolerance around limit value  m<sup>3</sup>/h

Reaction rate of diversity control  [0.01 ... 100] % per s



## Room exhaust air optimisation

Sum of technical V<sub>min</sub> of all Room exhaust air controller  m<sup>3</sup>/h

Tolerance for switching off the Room exhaust air controller  %

Delay time for switching off the Room exhaust air controller  s

## External volume flow shift via LonWorks<sup>®</sup>/BACnet/Modbus

Adjusting room air exchange rate (e.g. for external temperature control)

Characteristic	min	0 %	equals	min	<input type="text" value="0"/>	m <sup>3</sup> /h	(Minimum volume flow shift)
	max	100 %	equals	max	<input type="text" value="0"/>	m <sup>3</sup> /h	(Maximum volume flow shift)

# Room Controller config

## Pressure control - Pressure measuring point (AI5)

Characteristic min	<input type="text" value="0.00"/>	Volt	equals	min	<input type="text" value="-50"/>	Pa	Signal smoothing (1-50)	<input type="text" value="50"/>
max	<input type="text" value="10.00"/>	Volt	equals	max	<input type="text" value="50"/>	Pa		

## Pressure control - Parameters of pressure setpoint 1

Pressure setpoint	<input type="text" value="30"/>	Pa	Control tolerance	<input type="text" value="5"/>	Pa
Supply-Exhaust air balance	<input type="text" value="100"/>	m <sup>3</sup> /h			

(Positive value for overpressure / Negative value for depression)

Automatic calculation of limit value to achieve maximum actuating value

Increase pressure	<input type="text" value="15"/>	Pa
Decrease pressure	<input type="text" value="15"/>	Pa

Manual adaption of limit to achieve maximum actuating value.

Increase pressure

Minimum actuating value	<input type="text" value="0.30"/>	[0.01 ... 100]	Maximum actuating value	<input type="text" value="7.00"/>	[0.01 ... 100]
-------------------------	-----------------------------------	----------------	-------------------------	-----------------------------------	----------------

Decrease pressure

Minimum actuating value	<input type="text" value="0.60"/>	[0.01 ... 100]	Maximum actuating value	<input type="text" value="14.00"/>	[0.01 ... 100]
-------------------------	-----------------------------------	----------------	-------------------------	------------------------------------	----------------

Shift limits	min	<input type="text" value="0"/>	m <sup>3</sup> /h	Recommended values	-50	m <sup>3</sup> /h
	max	<input type="text" value="0"/>	m <sup>3</sup> /h		50	m <sup>3</sup> /h

## Alarm settings RMF alarm

- F1 Room alarms based on category 1 (Smoke extraction / UPS Battery operation)  Activate monitoring
- F2 Room alarms based on category 2 (Volume flow control)  Activate monitoring (Standard mode, High mode)  
 Activate monitoring (Low mode)
- F3 Room alarms based on category 3 (Hardware failure)  Activate monitoring

A7 Minimum Total exhaust air too low

Activate monitoring

Minimum Total exhaust air


Tolerance

%

Alarm delay time

s





Device	I/O	Volume flows	Operation modes	System	Alarm	Loupe	Damper positions	Terminal assignment	
<b>Within system detected EASYLAB devices</b>									
Fume cupboards				1				No. of devices with Room-Management-Function (RMF)	1
Room exhaust controller				1				Pressure control configured	no
Room supply controller				1				Number within system detected CL Terminations	2
TROX Adaption Module (TAM)				0				CL-Termination on this controller	active
Room supply controller for pressure control				0				System conflict	no
Room exhaust controller for pressure control				0					
Total number				3					



# Room Controller Diagnostics

Device I/O **Volume flows** Operation modes System Alarm Loupe Damper positions Terminal assignment

### Detailed information about integrated volume flows on this fume cupboard controller


Integrated volume flows	Exhaust air	Supply air	
DI 1 not configured as constant volume flow			
DI 2 not configured as constant volume flow			
DI 3 not configured as constant volume flow			
DI 4 not configured as constant volume flow			
DI 5 not configured as constant volume flow			
DI 6 not configured as constant volume flow			
AI 2 not configured as variable volume flow			
AI 3 not configured as variable volume flow			
AI 4 not configured as variable volume flow			
AI 5 not configured as variable volume flow			
LonWorks®/BACnet/Modbus	0	0	m³/h
Sum	0	0	m³/h

### Fume cupboard controller

Setpoint volume flow (Diversity factor = 100%)	432	m³/h
Setpoint volume flow (Diversity factor <= 100%)	432	m³/h
Current volume flow	414	m³/h
Damper position	39	% opened

### System

Setpoint Total exhaust air	400	m³/h
Current Total exhaust air	652	m³/h
Setpoint Total supply air	652	m³/h
Current Total supply air	670	m³/h
Total integrated exhaust air	0	m³/h
Total integrated supply air	0	m³/h
Total volume flow Fume cupboards	410	m³/h
Total volume flow Room exhaust controller	241	m³/h
Total volume flow Room supply controller	670	m³/h
Constant exhaust air	0	m³/h
Constant supply air	0	m³/h
Configured balance (Supply - Exhaust air)	0	m³/h
Volume flow shift - Pressure	0	m³/h
Volume flow shift - Temperature	0	m³/h



Device	I/O	Volume flows	Operation modes	System	Alarm	Loupe	Damper positions	Terminal assignment
--------	-----	--------------	-----------------	--------	-------	-------	------------------	---------------------

**Detailed information about controller damper positions in system**



	Damper position min [%]	Controller in operation mode Shut off	Damper position max [%]	Controller in operation mode Open
Fume cupboards	31	not present	31	not present
Room exhaust controller	35	not present	35	not present
Fume cupboards/Room exhaust controller	31	not present	35	not present
Room supply controller	54	not present	54	not present



### Room-Example 1

1. **Fume hoods with variable volume flow**
2. **Cupboards with constant 24h exhaust air**
3. **Room supply air is following the total exhaust**
4. **Room not airtight: No sealings on the doors**
5. **Room underpressure by volume flow difference**


## Room-Example 1



### Room-Example 2

1. Fume hoods with variable volume flow
2. Room exhaust is working vice versa to the hoods
3. Cupboards with constant 24h exhaust air
4. Room supply air is following the total exhaust
5. Room not airtight: No sealings on the doors
6. Room underpressure by volume flow difference

## Room-Example 2




### Room-Example 3

1. Fume hoods with variable volume flow
2. Room exhaust is working vice versa to the hoods
3. **Local exhaust with ON / OFF**
4. **functionSnorkels with manual dampers**
5. Cupboards with constant 24h exhaust air
6. Room supply air is following the total exhaust
7. Room not airtight: No sealings on the doors
8. Room underpressure by volume flow difference





## Room-Example 3



### Room-Example 4


1. Many fume hoods with variable volume flow
2. Room exhaust is working vice versa to the hoods
3. Local exhaust with ON / OFF function
4. Snorkels with manual dampers
5. Cupboards with constant 24h exhaust air
6. **Two Room supply air is following the total exhaust (devided air 50% / 50%)**
7. Room more airtight: Sealings on the doors
8. Room is pressure controlled  
volume flow / pressure cascade

# EASYLAB – room example



**Room-Example 4**


**With EASYLAB**  
**automatically possible**



### Room-Example 5

1. Fume hoods with variable volume flow
2. Room exhaust is working vice versa to the hoods
3. Local exhaust with ON / OFF function
4. Snorkels with manual dampers
5. Cupboards with constant 24h exhaust air
6. Room supply air is following the total exhaust
7. Room more airtight: Sealings on the doors
8. Room is pressure controlled  
volume flow / pressure cascade

Not so easy to keep him on the same size



## Physical Background

Roomleakage: 0,001 m<sup>2</sup>

About 1 mm gap below the door

Room pressure = 463 Pa  
100 m<sup>3</sup>/h Vol.-Difference

**Guilty!**



According to Bernoulli:

$$\Delta p = \frac{\rho}{2} * \left( \frac{\dot{V}}{A} \right)^2$$

$\dot{V}$  = Volume flow difference

A = Roomleakage

$\rho$  density of air 1,2 kg/m<sup>3</sup>


If the area is nearly zero this part of the formular is nearly unlimited!





Then little changes in the volume flow have extreme effects

The volume flow difference is independent of the room size!

about 1mm gap

$$\Delta p = \frac{\rho}{2} * \left( \frac{\dot{V}}{A * \mu} \right)^2$$







## Room-Example 5

